



COBALT HEALTH BARIATRIC AIR PRESSURE CARE CHAIR



SKU: C01-BARIATRIC-BLU

PRODUCT DESCRIPTION

Purpose & Benefits

The Cobalt range of Air Chairs are specifically designed for use in all environments where users require elevated levels of comfort, coupled with low - high risk pressure care and mobility.

The Bariatric Air Chair features a versatile multi positional design with a 60cm wide seat and an increased 200kg SWL. Each Air Chair employs a combination of soft high-density foam, adjustable air cells (seat & lumbar sections) and a multi-directional stretch cover.

Every chair is highly adjustable, with independently operated tilt-in-space, recline and leg raise functions that provide a wide array of repositioning and pressure redistribution options.

Features

- Highly adjustable multi positional pressure care seating system
- Available in Slimline, Classic and Bariatric sizes
- Independently operated Tilt-in-space, recline and leg raise functions
- Adjustable air cells in seat and lumbar sections
- Drop down arm rests & fold back wings assist with cleaning and transfers
- Adjustable & removable Head Rest
- "No Gap" seating surface
- 4 x large 150mm swivel castors to aid mobility
- 3 x Non marking brake castors & 1 x steering castor
- Fluid resistant, vapour permeable, 4-way stretch cover
- Cover can be removed for cleaning or replaced if damaged
- Fire retardant cover and foam

- Height adjustable flip up footrests included with each Air Chair
- Optional Meal Tray available
- SWL 200kg
- Warranty 3 Years Frame, 1 Year Gas Struts, foam and air cells

Specifications

Specification	Details
Overall Dimensions	130 x 80 x 125 cm
Seat Width (between arms)	60 cm
Seat Depth	56 cm
Seat Height	70 cm
Seat Recline Angle	97-136°
Seat Tilt Angle	0-17°
Leg Rest Angle	0-90°
Available for Hire	Yes
Weight	60 kg
SWL (Safe Working Load)	200 kg
Warranty	3 Years Frame, 1 Year Gas struts, foam and air cells

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au