




## **COBALT HEALTH CLASSIC AIR PRESSURE CARE CHAIR**



**SKU:** C01-CLASSIC-BLU

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The Cobalt Health Air Chair Classic is designed to provide superior comfort and pressure relief for patients who require extended seated care. With a focus on both function and support, it combines deep, pressure-reducing foam with high backrest design, lateral supports, and adjustable leg elevation, making it ideal for hospital, aged care, or home environments.

Fitted with large rear wheels and front castors for smooth manoeuvrability, the chair can be positioned or reclined to suit individual needs and safely locked in place. It's an ideal choice for patients with low mobility, offering improved positioning, circulation support, and enhanced comfort for longer durations of sitting.

### **Features**

- High-back design with lateral support wings
- Deep pressure-relieving foam seat and backrest
- Tilt-in-space recline function
- Height-adjustable leg rest
- Padded armrests for added comfort

- Swing-away and removable leg supports
- Easy-rolling rear wheels and front castors
- Lockable castors for secure positioning
- Suitable for hospitals, aged care, or in-home use
- Designed for users with low mobility or postural care needs

## Specifications

- Safe Working Load: 150 kg
- Seat Height: 700 mm
- Seat Width: 510 mm
- Seat Depth: 560 mm
- Overall Height: 1300 mm
- Overall Width: 710 mm
- Overall Depth: 1250 mm
- Seat Back Adjustment Angle: 97° – 136°
- Footrest Adjustment Angle: 0° – 90°
- Seat Tilt Adjustment Angle: 0° – 17°
- Included Accessories: Headrest, Footrest Pads

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)