



## **COBALT HEALTH SLIMLINE AIR PRESSURE CARE CHAIR**



**SKU:** C01-SLIMLINE-BLU

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The Cobalt Health Air Chair Slimline offers the same therapeutic support and pressure care benefits as the Classic model, but in a narrower frame—perfect for smaller users or tighter clinical spaces. Designed to maximise comfort and safety, this chair includes tilt-in-space functionality, reclining backrest, and pressure-relieving foam to support long-term seated care.

With smooth-rolling castors, lockable wheels, and a lightweight frame, it's ideal for both residential and professional care environments. Adjustable components allow for tailored positioning, while included headrest and footrest pads offer additional comfort and support during extended use.

### **Features**

- Slimline frame for compact spaces and smaller users
- Tilt-in-space and recline functionality
- Padded lateral supports and high backrest
- Pressure-relieving seat and backrest foam
- Height-adjustable footrest with removable pads

- Lockable castors for safe positioning
- Lightweight and easy to manoeuvre
- Ideal for hospitals, aged care, or home use
- Includes headrest and footrest pads
- Optional meal tray available

## **Specifications**

- Safe Working Load: 150 kg
- Seat Height: 700 mm
- Seat Width: 450 mm
- Seat Depth: 560 mm
- Overall Height: 1300 mm
- Overall Width: 650 mm
- Overall Depth: 1250 mm
- Seat Back Adjustment Angle: 97° – 136°
- Footrest Adjustment Angle: 0° – 90°
- Seat Tilt Adjustment Angle: 0° – 17°
- Included Accessories: Headrest, Footrest Pads
- Optional Accessory: Meal Tray

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)