

## COBALT HEALTH BEDSIDE CABINET SINGLE LOCK



**SKU:** A02-CBT-BCH-SL

## PRODUCT DESCRIPTION

### Purpose and Benefits

The Bedside Cabinet with Single Lock Drawer in Beech finish is a practical and secure storage solution designed to complement the CS8 bed in both clinical and home care settings. It provides convenient access to personal items while keeping valuables and medications safely locked away.

Its compact design includes smooth-rolling concealed castors for easy repositioning, and its neutral finish blends well with modern care room environments. Ideal for aged care, hospitals, and assisted living setups, this cabinet is designed for everyday use around the bed.

- Keeps essentials within reach of the bed
- Lockable drawer secures medications and valuables
- Easy to move with discreet castors
- Matches the style and height of the CS8 bed and similar models

### Features

- Three-drawer configuration: two standard drawers plus one deep dumper drawer
- Top drawer includes integrated lock for secure storage

- Optional drawer dividers available for internal organisation
- Powder-coated steel body with chrome-plated handles
- Smooth-glide concealed castors for mobility
- Durable, wipe-clean surfaces suitable for medical environments
- Available in Beech and Cream finishes

## **Specifications**

- Dimensions: 45 cm (Length) x 45 cm (Width) x 90 cm (Height)
- Frame Material: Powder-coated steel
- Drawer Handles: Chrome-plated
- Finish: Beech
- Castors: 4 concealed wheels for discreet movement
- Locking: Single lock on top drawer
- Compatibility: Designed for use with CS8 and similar bed models