



COLOPLAST CATHETER VALVE

\$14.00



SKU: 13290100

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Intended to be connected to an indwelling catheter and is used to drain the urine either directly into the toilet/container or into a bag connected to the valve.

FEATURES

- Provides a discreet form of continence management.
- Helps maintain bladder tone and capacity by mimicking normal bladder function via periodic opening to drain the bladder.
- The valve can be used with or without a drainage bag.
- Ridged connector for secure connection and safe removal.
- The lever-action tap is easy to use.
- Suitable for use with all types of internal or external catheters.