

**COLOPLAST SPEEDICATH MALE CATHETER****\$9.00****SKU: 28414****PRODUCT DESCRIPTION****PURPOSE & BENEFITS**

SpeediCath catheters are safe, hygienic catheters that offer ease of use, greater comfort and minimal risk of urethral damage for all catheter users. Instantly ready to use without the need to add neither water or waiting for the coating to get activated, easy to open and simple to handle SpeediCath represents a more convenient way to catheterise. The uniform hydrophilic surface coating ensures safe catheterisation time after time.

FEATURES

- A smooth even coating minimises the risk of urethral damage
- The SpeediCath's are pre-hydrated and thus ready to use
- Intuitive and easy to use

SPECIFICATIONS

- Length - 40cm
- Profile - Straight Tip
- Size - FG10
- Gender - Male