




HEMOCRAFT ALUMINIUM DELUXE BEDSIDE COMMODE

\$347.00



SKU: NOV-COM71137

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Homecraft Aluminium Deluxe Bedside Commode is comfortable and hygienic yet modest. Features rust resistant aluminium frame, which is sturdy and height adjustable. This commode chair features non-slip rubber tips and an easily removable backrest. The lift out bowl with lid, padded vinyl contoured backrest, toilet seat and cover offer extra comfort and pressure relief whilst maintaining ease of cleaning and hygiene.

FEATURES

- Non-rustable aluminium frame
- Height adjustable
- Non-slip rubber tips
- Removable backrest
- Lift out bowl with lid
- Padded vinyl contoured backrest, toilet seat and cover

SPECIFICATIONS

- Seat Width Between Arms: 450mm
- Overall Width & Depth: 550mm
- Seat Height Range: 440mm - 590mm
- Armrest Height (Seat to Top of Handle): 180mm
- Maximum User Weight: 150kg

