



## COMSEAT BARIATRIC ALBANY HIGH BACK CHAIR



**SKU:** 400246

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Albany High Back Chair** is a durable, Australian-made seating solution designed for comfort, support and longevity in commercial and healthcare environments. Engineered from heavy-duty Australian steel and reinforced timber, it provides exceptional strength and stability for users up to 300 kg. Its high back design, deep cushioning and modern high-resilience foam deliver superior comfort, making it ideal for long-term use in aged care, hospitals, waiting areas and home settings.

### Features

- High back design provides superior postural support and comfort
- Heavy-duty frame constructed from broad diameter Australian steel tubing
- Steel-braced timber seat and back for enhanced durability and stability
- Independently certified to a Safe Working Load (SWL) of 300 kg
- Deep-cushioned seat with Australian-made high resilience foam for long-lasting comfort
- Available in a wide range of colours and fabrics to suit different interiors
- Optional matching foot stool available
- Proudly made in Australia to high commercial-grade standards

## Specifications

Specification	Details
Model	Albany High Back Chair
Safe Working Load (SWL)	300 kg (independently certified)
Frame Material	Heavy-duty Australian steel tubing
Seat & Back	Steel-braced timber construction
Cushioning	High resilience Australian foam for deep cushioned comfort
Customisation	Available in various colours and fabrics
Accessories	Optional matching foot stool
Country of Manufacture	Australia

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriiffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
 1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)

