



#### **COMSEAT BARIATRIC ALBANY HIGH BACK CHAIR**



**5KU:** 400246

#### PRODUCT DESCRIPTION

## **Purpose & Benefits**

The **Albany High Back Chair** is a durable, Australian-made seating solution designed for comfort, support and longevity in commercial and healthcare environments. Engineered from heavy-duty Australian steel and reinforced timber, it provides exceptional strength and stability for users up to 300 kg. Its high back design, deep cushioning and modern high-resilience foam deliver superior comfort, making it ideal for long-term use in aged care, hospitals, waiting areas and home settings.

### **Features**

- High back design provides superior postural support and comfort
- Heavy-duty frame constructed from broad diameter Australian steel tubing
- Steel-braced timber seat and back for enhanced durability and stability
- Independently certified to a Safe Working Load (SWL) of 300 kg
- Deep-cushioned seat with Australian-made high resilience foam for long-lasting comfort
- Available in a wide range of colours and fabrics to suit different interiors
- Optional matching foot stool available
- Proudly made in Australia to high commercial-grade standards



# **Specifications**

**Specification** Details

Model Albany High Back Chair

Safe Working Load (SWL) 300 kg (independently certified)
Frame Material Heavy-duty Australian steel tubing

Seat & Back Steel-braced timber construction

Cushioning High resilience Australian foam for deep cushioned comfort

Customisation Available in various colours and fabrics

Accessories Optional matching foot stool

Country of Manufacture Australia