



CONFIGURA ADVANCE MANUAL CARE CHAIR



SKU: CR5438

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Configura Advance Manual Care Chair, a highly configurable specialist chair designed to meet complex and challenging seating needs without breaking the bank. This innovative and cost-effective seating solution combines state-of-the-art features and standards that set it apart from other chairs in its class.

FEATURES

- **37 Degrees of Tilt in Space & Full-Length Compensating Legrest:** Experience unparalleled comfort and flexibility with a wide range of tilt options, along with a compensating legrest.
- **Lockable, Reduced Shear Backrest Recline:** Maximize pelvic stability with a backrest recline feature that minimizes shear forces, giving you the support you need.
- **Tool-less Adjustments with ClickLock Adjustment System:** Say goodbye to cumbersome tools and enjoy easy and hassle-free adjustments with the innovative ClickLock system.
- **Upgradeable Backrest Options:** Tailor your chair to your specific needs with customizable backrest options.
- **Removable, Angle, and Height-Adjustable Footplate as Standard:** For personalized comfort and fit, this chair comes with a footplate that can be adjusted in height and angle or entirely removed.
- **Visco Integrated Self-Supply Cushion as Standard:** Measuring 17.5 x 17.5 x 3 inches, the self-supply cushion ensures additional comfort and relaxation.
- **Drop-Down Armrest Option:** Facilitating easier transfers and sling insertions, the drop-down armrest is a functional feature that adds to the chair's versatility.

- **Forward Tilt Option (Standard on Electrically Operated Chairs):** Benefit from the added forward tilt feature for a more ergonomic experience.
- **Electric Option Available:** For those seeking automation and ease of use, the Configura Advance is also available in an electric version.

SPECIFICATIONS

Specification	Details
Colour	Charcoal
Fabric	Vinyl fabric with VP on all contact areas
Seat Width Options	350 mm (14 in) – 550 mm (22 in)
Seat Depth Options	400 mm (16 in) – 550 mm (22 in)
Seat Height Options	330 mm (13 in) – 530 mm (21 in)
Tilt-in-Space	37 degrees
Backrest Recline	130 degrees
Forward Tilt	10 degrees
Weight Capacity	160 kg (25 st)
Mechanism	Manual tilt-in-space

VIDEO

[embed]<https://www.youtube.com/watch?v=msXT6VjWgoU>[/embed]

BROCHURE

[Configura Advance Brochure](#)