



## CONFIGURA ADVANCE MANUAL CARE CHAIR



**SKU:** CR5438

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Configura Advance Manual Care Chair, a highly configurable specialist chair designed to meet complex and challenging seating needs without breaking the bank. This innovative and cost-effective seating solution combines state-of-the-art features and standards that set it apart from other chairs in its class.

### FEATURES

- **37 Degrees of Tilt in Space & Full-Length Compensating Legrest:** Experience unparalleled comfort and flexibility with a wide range of tilt options, along with a compensating legrest.
- **Lockable, Reduced Shear Backrest Recline:** Maximize pelvic stability with a backrest recline feature that minimizes shear forces, giving you the support you need.
- **Tool-less Adjustments with ClickLock Adjustment System:** Say goodbye to cumbersome tools and enjoy easy and hassle-free adjustments with the innovative ClickLock system.
- **Upgradeable Backrest Options:** Tailor your chair to your specific needs with customizable backrest options.
- **Removable, Angle, and Height-Adjustable Footplate as Standard:** For personalized comfort and fit, this chair comes with a footplate that can be adjusted in height and angle or entirely removed.
- **Visco Integrated Self-Supply Cushion as Standard:** Measuring 17.5 x 17.5 x 3 inches, the self-supply cushion ensures additional comfort and relaxation.
- **Drop-Down Armrest Option:** Facilitating easier transfers and sling insertions, the drop-down armrest is a functional feature that adds to the chair's versatility.

- **Forward Tilt Option (Standard on Electrically Operated Chairs):** Benefit from the added forward tilt feature for a more ergonomic experience.
- **Electric Option Available:** For those seeking automation and ease of use, the Configura Advance is also available in an electric version.

## SPECIFICATIONS

Specification	Details
Colour	Charcoal
Fabric	Vinyl fabric with VP on all contact areas
Seat Width Options	350 mm (14 in) – 550 mm (22 in)
Seat Depth Options	400 mm (16 in) – 550 mm (22 in)
Seat Height Options	330 mm (13 in) – 530 mm (21 in)
Tilt-in-Space	37 degrees
Backrest Recline	130 degrees
Forward Tilt	10 degrees
Weight Capacity	160 kg (25 st)
Mechanism	Manual tilt-in-space

## VIDEO

[embed]https://www.youtube.com/watch?v=msXT6VjWgoU[/embed]

## BROCHURE

[Configura Advance Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)