

**CONFIGURA BARIATRIC CHAIR, BLACK VINYL****\$11,963.26****SKU: CR5472****PRODUCT DESCRIPTION****Purpose & Benefits**

The **Configura Bariatric Chair** is a high-capacity electric riser recliner designed specifically for bariatric users. It features Configura's trusted adjustability and support technology in a durable, clinical-grade format. This model includes independent motors for the leg-rest, backrest, and tilt-in-space (TIS) functions, allowing precise postural control and improved comfort over long periods of use.

Upholstered in black vinyl for easy cleaning and infection control, the chair also includes the 28-inch arm width kit, making it suitable for broader users requiring greater seat space and stability.

Ideal for residential care, hospitals, and home use, it supports users with reduced mobility and those who need enhanced pressure management and postural assistance.

**Features**

- Bariatric-rated electric riser recliner chair
- Black vinyl upholstery for hygiene and durability
- Independent motors: leg-rest, backrest, and TIS
- High weight capacity for safety and support

- Supplied with 28" arm kit
- Adjustable seat depth (20" or 22")
- Configura modular technology for customisation

## Specifications

- **Weight Capacity:** 254kg
- **Leg-rest Weight Capacity:** 127kg
- **Seat Height:** 16" (406mm)
- **Seat Width:** 28" (711mm)