



CONFIGURA BARIATRIC CHAIR - BLACK VINYL



SKU: CR5472-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **Configura Bariatric Chair** is a high-capacity electric riser recliner designed specifically for bariatric users. It features Configura's trusted adjustability and support technology in a durable, clinical-grade format. This model includes independent motors for the leg-rest, backrest, and tilt-in-space (TIS) functions, allowing precise postural control and improved comfort over long periods of use.

Upholstered in black vinyl for easy cleaning and infection control, the chair also includes the 28-inch arm width kit, making it suitable for broader users requiring greater seat space and stability.

Ideal for residential care, hospitals, and home use, it supports users with reduced mobility and those who need enhanced pressure management and postural assistance.

Features

- Bariatric-rated electric riser recliner chair
- Black vinyl upholstery for hygiene and durability
- Independent motors: leg-rest, backrest, and TIS
- · High weight capacity for safety and support



- Supplied with 28" arm kit
- Adjustable seat depth (20" or 22")
- Configura modular technology for customisation

Specifications

• Weight Capacity: 254kg

• Leg-rest Weight Capacity: 127kg

• **Seat Height**: 16" (406mm)

• Seat Width: 28" (711mm)

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au



NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au