




CONFIGURA BARIATRIC RISER RECLINER LOUNGE CHAIR

\$8,666.00 - \$10,358.00



SKU: CR5412-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Configura Bariatric Riser Recliner Chair: Unmatched Adaptability and Customization in One Revolutionary Lift Chair!

The Configura Bariatric Riser Recliner Chair is a game-changer in the world of lift chairs, offering unrivaled adaptability and a wide range of custom options. Designed with user comfort and satisfaction in mind, this chair is a true innovation.

Experience the Configura Bariatric Riser Recliner Chair in Duratek fabric, a durable and high-quality material. With a Tilt-in-Space of 16", this chair provides exceptional support and positioning. The chair comes equipped with a Visco cushion and a Multi-Adjustable Pillow Backrest, ensuring optimum comfort and postural support.

Our chairs are meticulously designed to provide comfortable seating for users with mild to moderate postural requirements. They are specifically crafted to help reduce the formation of pressure ulcers and assist users with limited mobility in standing or sitting.

What sets the Configura Bariatric Riser Recliner Chair apart is its incredible adjustability. With three different sizes available, we guarantee the perfect fit for each user. Our chairs are thoughtfully designed to ensure maximum comfort and support, tailored to individual needs.

[CLICK HERE TO VIEW THE FULL BROCHURE](#)

Configura Bariatric is appropriate for:

- Bariatric users
- Users at risk of pressure damage
- Users with current pressure damage (subject to professional assessment)
- Ambulant, semi-ambulant and non-ambulant users
- Accommodation and encouragement of good posture
- Users with Neurological conditions
- Users with a kyphosis
- Users with oedematous legs and/or fluid build-up
- The flat or vertical rise function is particularly suited to users with Muscular Dystrophy who need to lock knees before transferring

It is not appropriate for:

- Demanding postural requirements
- Users with fixtures and contractures
- Users with a limited extension at the knees
- Users who need to be transported from room to room

FEATURES

- The Dual motor design allows independent adjustment of leg rest elevation and backrest recline.
- Chair width can be adjusted easily via the removable armrests, allowing three different width options in the same chair.
- Breathable and easily wipeable fabric on seat base and backrest
- The adjustable nature of the chair means that even if circumstances change and pressure care risk increases or other adaptations are required, these needs can often be met without having to replace the chair

SPECIFICATIONS

Minimum Overall Height (mm)	1150
Maximum Overall Height (mm)	1150
Cover Material	Duratek
Colour	Sandstone
Warranty (months)	24
DVA Approved (if applicable)	Yes
Safety Accreditation (if applicable)	BS EN 12520:2010 EN 55014
Maximum Seat Height (mm)	406
Seat Width (mm)	609 711
Seat Depth (mm)	508 550
Backrest Height (from Seat) (mm)	700

Number of Motors	Dual Motor
Length Fully Reclined (mm)	1700
Bulky Product	No
Assembly Required	No
Available For Hire	No
Quote Only	No