

CONFIGURA BARIATRIC PROFILED HEADREST





5KU: CA2448

PRODUCT DESCRIPTION

Purpose & Benefits

The **Configura Bariatric Profiled Headrest (VP Material)** provides increased head support for users with reduced head control. It forms the top cushion of the Waterfall Backrest and can be removed or replaced as needed.

Features

- Top cushion of the Waterfall Backrest
- Removable and replaceable design
- Provides enhanced head support
- Suitable for users with reduced head control
- Compatible with all Configura chairs

Specifications

• Fabric: VP Material

• Compatibility: All Configura chairs

