




## **CONFIGURA COMFORT DURATEK RISER RECLINER LOUNGE CHAIR**



**SKU:** CR540-H

### **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Configura Comfort is a revolutionary new chair with unparalleled adaptability and custom options all in one chair!

The chairs are designed to provide comfortable seating to users with mild to moderate postural requirements, to help reduce the formation of pressure ulcers and help to stand or sit users with limited mobility. The higher position of the leg rest when raised is ideal for assisting with swelling and circulation issues in the legs.

The chair offers great adjustability, in 3 different sizes, to ensure the chair fits the user perfectly.

In addition, there are 3 integrated options for the lift function:

1. Standard rise lift with 23° of tilt
2. Medium rise lift with 14° of tilt
3. Flat (or Vertical) rise with 2° of tilt and 10" (260mm) of rise (this function is not available on the small model)

### **FEATURES**

- Designed to provide comfortable seating to users with mild to moderate postural requirements, to help reduce the formation of pressure ulcers
- 'Reduced shear' to maximise pelvic stability
- Pressure-reducing fabric on the backrest and armrest as standard.
- Castors for easy manoeuvrability

- Integrated self-supply overlay for easy integration of existing cushions

The Configura Comfort is ideal for:

- Users at risk of pressure damage
- Users with current pressure damage (subject to professional assessment)
- Ambulant, semi-ambulant and non-ambulant users
- Accommodation and encouragement of good posture
- Users with Neurological conditions
- Users with a kyphosis
- Users with oedematous legs and/or fluid build-up
- The flat or vertical rise function is particularly suited to users with Muscular Dystrophy who need to lock knees before transferring

## SPECIFICATIONS

Chair Size	Part Number	Seat Depth	Seat Height	Seat Width	Weight Capacity
SMALL	CR5404	457 mm (18"), 431 mm (17")*	406 mm (16"), 508 mm (20")	457 mm (18"), 457 mm (18")**	130 kg
MEDIUM	CR5406	457 mm (18"), 482 mm (19")*	457 mm (18")	508 mm (20"), 508 mm (20")**	160 kg
LARGE	CR5406	457 mm (18")	457 mm (18")	559 mm (22")	160 kg

## BROCHURE

[Click here to view full brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)