




## **CONFIGURA COMFORT BLACK VINYL RISER RECLINER LOUNGE CHAIR**

**\$4,894.00 - \$5,499.00** Price range: \$4,894.00 through \$5,499.00



**SKU:** CR5444-H

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Configura Comfort electric recliner lift chair is designed to provide comfortable seating to users with mild to moderate postural requirements, to help reduce the formation of pressure ulcers and help to stand or sit users with limited mobility.

All Configura chairs have tilt-in-space mechanisms, integrated, interchangeable pressure care and interchangeable posture management options. If you require adjustments or alterations to the basic size models, please contact or visit us at your nearest store.

### **Appropriate for:**

- Users at risk of pressure damage
- Users requiring assistance sitting and rising from a seated position safely and comfortably
- Encouraging good posture

### **FEATURES**

- 'Reduced shear' recline to maximise pelvic stability
- Pressure reducing fabric on backrest and armrest as standard
- Castors for easy manoeuvrability
- Integrated self-supply overlay for easy integration of existing cushion

**[CLICK HERE TO VIEW THE FULL BROCHURE](#)**

## **SPECIFICATIONS**

- Small
  - Floor to seat height – 40.6 cm
  - Seat width – 45.7 cm
  - Seat depth – 45.7 cm
- Medium
  - Floor to seat height – 45.7 cm
  - Seat width – 50.8 cm
  - Seat depth – 45.7 cm
- Large
  - Floor to seat height – 45.7 cm
  - Seat width – 55.9 cm
  - Seat depth – 45.7 cm

*\*dimension adjustments available on request.*