



CONFIGURA COMFORT VINYL RISER RECLINER LOUNGE CHAIR



SKU: CR5426-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Configura Comfort electric recliner lift chair is designed to provide comfortable seating to users with mild to moderate postural requirements, to help reduce the formation of pressure ulcers and help to stand or sit users with limited mobility.

All Configura chairs have tilt-in-space mechanisms, integrated, interchangeable pressure care and interchangeable posture management options. If you require adjustments or alterations to the basic size models, please contact or visit us at your nearest store.

Appropriate for:

- Users at risk of pressure damage
- Users requiring assistance sitting and rising from a seated position safely and comfortably
- Encouraging good posture

FEATURES

- 'Reduced shear' recline to maximise pelvic stability
- Pressure reducing fabric on backrest and armrest as standard
- Castors for easy manoeuvrability
- Integrated self-supply overlay for easy integration of existing cushion

[CLICK HERE TO VIEW THE FULL BROCHURE](#)

SPECIFICATIONS

- Small
 - Floor to seat height - 40.6 cm
 - Seat width - 45.7 cm
 - Seat depth - 45.7 cm
 - SWL: 160kg
- Medium
 - Floor to seat height - 45.7 cm
 - Seat width - 50.8 cm
 - Seat depth - 45.7 cm
 - SWL: 160kg
- Large
 - Floor to seat height - 45.7 cm
 - Seat width - 55.9 cm
 - Seat depth - 45.7 cm
 - SWL: 160kg

**dimension adjustments available on request.*

**Please note: May require assembly.*

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au