



ERGONOMIC FOREARM CRUTCHES

\$38.00 Original price was: \$38.00.<u>\$34.00</u>Current price is: \$34.00.



SKU: JAN-124-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

These Ergonomic Forearm Crutches are Ergonomically shaped anatomical handgrip forearm crutch designed to distribute pressure more evenly across the hand.

KEY FEATURES

- Durable and strong, these Ergonomic Grip Forearm Crutches are constructed from high tensile extruded aluminium making them both lightweight and very strong.
- Ergonomic palm grip offers pressure distribution for maximum comfort.
- Fitted with high quality rubber tips for stability.
- Both hand to floor and palm grip to arm cuff sections are adjustable.
- Maximum User Weight 160kg

SPECIFICATIONS

- M (Hand Grip Height 66 94cm, Arm Ring Adjustment- 20.5 -27.5cm)
- L (Hand Grip Height 83 111cm, Arm Ring Adjustment- 20.5 -27.5cm)

Buy the Ergonomic Forearm Crutches Australia online at MediMart using Afterpay or NDIS method today!