



CURA1 CORDLESS TRANSMITTER/BATTERY PACK FOR WIRELESS SENSOR PADS & MATS



SKU: 2719

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Cura1 Wireless Pad Transmitter connects with the Cura1 [Floor](#), [Bed](#) and [Chair](#) pads. You will also need the [Cordless Falls Monitor](#).

When a resident is placed on a bed or chair mat the pad transmitter will detect their weight. Therefore any weight removed from the pad will trigger an alert signal via A-CBM and as a result, send an out of bed alarm to notify a carer.

The transmitter uses a non- replaceable long life battery (Lith-ION) with 12 month life expectancy covering approximately 20-50 m range.

This new model now incorporates “sleep mode” when not in use therefore extending battery life. Battery condition is monitored via the bedside control unit (Cura1 Wireless Monitor).

Up to 5 transmitters can be in use at the one time with one monitor and as a result can have a bed, chair and floor mat in operation simultaneously.

All Cura1 pads and floor mats have the the same unique connection and it is for that reason that the transmitter can be swapped between pads and mats as required.

Above all due to the size of the transmitter is not easily turned off by the elderly or aged due to the dexterity required.

FEATURES

- 2mm x1.2mm 4pin pitch plug

- 433MHz operating frequency
- Between 20- 50m range*
- Uses a Lith-ION long life battery – 1 year battery life 6 month warranty
- Will go into sleep mode when not in use to save battery life
- Status/in-use LED indication
- Communicates with A- CBM- bedside controller
- Incorporates on/off
- Fits within floor, bed & chair sensor transmitter enclosure

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au