




CURA1 2700 SIX CHANNEL CORDLESS FALLS MONITOR



SKU: 2700

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Cura1 Wireless bed pad kit home use is an industry standard falls prevention system, that allows for a carer in the home to have more flexibility and be alerted if an unassisted bed exit is detected.

Simple operation and maintenance.

No more cords between monitor and pads.

Can be used with up to six cordless sensor pads or other triggering devices.

PRODUCT HELP

To use the Cura Falls Monitor to its full potential, you will need to also purchase the [Cordless Transmitter](#), and either or all of the [Bed Exit Pad](#), [Chair Exit Pad](#), and the [Floor Exit Pad](#). All these products are available for individual purchase or you can simply buy them in a bundle from this page.

FEATURES

- Modern ergonomic design
- Monitor up to 7 devices simultaneously - a single Premium pad/mat and up to 6 wireless devices.
- Comes with durable ABS plastic enclosure and Caregiver Reset Key.
- Alarm memory and auto re-arm.

- Amber LED visual alarm indicator.
- Caregiver secure reset option.
- Battery operated or optional external plug pack.
- Up to 6 months battery life (nominal).
- Easy to access battery compartment.
- Can be connected to a nurse call system.
- 3 selectable alarm tones.
- Easy thumbwheel adjustable alarm volume.
- “Pause” feature allows caregiver to attend to the patient without setting the alarm off.
- Alerts for lost nurse call connection, lost AC power supply & low battery.

Safety Features

- Non resettable “Pad Signal Lost” warning.
- Caregiver secure reset option.
- All safety alerts sound even if alarm volume set to silent, including low battery, AC lost, pad signal lost, nurse call lost & change pad.
- Cable strain relief & breakaway cord technology prevents damage & helps reduce risk of entanglement when used with Premium Pads & Mats.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au