




CURA1 CORDLESS BED PAD



SKU: CUR-2629

PRODUCT DESCRIPTION

Purpose & Benefits of the Cura1 Cordless Bed Pad

The Cura1 Cordless Bed Pad is designed to enhance patient safety by alerting caregivers when an individual attempts to leave their bed. Its cordless design minimizes tripping hazards and simplifies setup, making it an essential tool in various care settings.

Features

- **Cordless Operation:** Eliminates the need for cords between the bed pad and monitor, reducing tripping hazards and allowing for flexible placement.
- **Easy to Clean:** Made from non-allergenic soft PVC, the pad is simple to wipe down, ensuring hygiene standards are maintained.
- **Moisture Resistance:** Equipped with vent holes in the neck area to reduce the risk of moisture ingress, enhancing durability.
- **Reliable Performance:** Designed to provide consistent monitoring with minimal false alarms over its lifespan.

Specifications

- **Dimensions:** 760 mm x 510 mm
- **Material:** Non-allergenic soft PVC
- **Weight Capacity Range:** 30 kg to 200 kg for reliable operation
- **Compatibility:** Requires the Cura1 2719 Cordless Pad Transmitter for operation
- **Warranty:** Guaranteed 1-year warranty

The Cura1 Cordless Bed Pad is a practical solution for caregivers seeking to monitor patient movement discreetly and effectively, thereby reducing the risk of falls and enhancing overall safety.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au