



CURA1 CORDLESS BED PAD

\$141.00



SKU: CUR-2629

PRODUCT DESCRIPTION

Purpose & Benefits of the Cura1 Cordless Bed Pad

The Cura1 Cordless Bed Pad is designed to enhance patient safety by alerting caregivers when an individual attempts to leave their bed. Its cordless design minimizes tripping hazards and simplifies setup, making it an essential tool in various care settings.

Features

- **Cordless Operation:** Eliminates the need for cords between the bed pad and monitor, reducing tripping hazards and allowing for flexible placement.
- **Easy to Clean:** Made from non-allergenic soft PVC, the pad is simple to wipe down, ensuring hygiene standards are maintained.
- **Moisture Resistance:** Equipped with vent holes in the neck area to reduce the risk of moisture ingress, enhancing durability.
- **Reliable Performance:** Designed to provide consistent monitoring with minimal false alarms over its lifespan.

Specifications

- **Dimensions:** 760 mm x 510 mm
- **Material:** Non-allergenic soft PVC
- **Weight Capacity Range:** 30 kg to 200 kg for reliable operation
- **Compatibility:** Requires the Cura1 2719 Cordless Pad Transmitter for operation
- **Warranty:** Guaranteed 1-year warranty

The Cura1 Cordless Bed Pad is a practical solution for caregivers seeking to monitor patient movement discreetly and effectively, thereby reducing the risk of falls and enhancing overall safety.