



CURA1 CORDLESS CHAIR PAD

\$112.00 Original price was: \$112.00.\$101.00 Current price is: \$101.00.



SKU: 2628

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

This wireless chair pad is used on chairs, recliners, wheel chairs etc. of those at the risk of having a fall.

When used in conjunction with a Cura1 Transmitter and a Cura1 wireless monitoring device, the removal of pressure generates an alarm.

You will need to have the <u>Cural Transmitter</u> and <u>Cural Cordless Falls Monitor</u> for everything to work.

Like the bed sensor it is an easy to clean non-allergenic soft PVC material printed with detailed instructions.

Please note: This is the chair pad ONLY, used for replacement parts and spares.

FEATURES

- Simple operation and maintenance
- Affordable price
- · Easy to clean
- For use with cural 2719 Cordless Pad Transmitter
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress