



## CURA1 CORDLESS CHAIR PAD

**\$112.00**



**SKU: 2628**

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

This wireless chair pad is used on chairs, recliners, wheel chairs etc. of those at the risk of having a fall.

When used in conjunction with a Cura1 Transmitter and a Cura1 wireless monitoring device, the removal of pressure generates an alarm.

You will need to have the [Cura1 Transmitter](#) and [Cura1 Cordless Falls Monitor](#) for everything to work.

Like the bed sensor it is an easy to clean non-allergenic soft PVC material printed with detailed instructions.

Please note: This is the chair pad ONLY, used for replacement parts and spares.

### FEATURES

- Simple operation and maintenance
- Affordable price
- Easy to clean
- For use with cura1 2719 Cordless Pad Transmitter
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress