


CURA1 2624 CORDLESS FLOOR MAT



SKU: 2624B

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Cura1 Wireless Floor Mat has the benefit of wireless out of bed detection methods in that there are no messy wires or potential trip hazards. As the resident places their foot on this floor mat the carer will be notified via the facilities nurse call system. This Cura1 brand mat has a “yoga” soft top and is waterproof and washable. It is low profile and made of a non-slip material. It has a press-stud protective pocket to house the transmitter/battery pack (sold separately).

Wireless floor mats and transmitter/battery packs are sold separately.

The transmitter/battery pack which is required for this mat is the [Cura1 Transmitter Module](#) and the Monitor is the [Cura1 Cordless Falls Monitor](#).

KEY FEATURES

- Simple operation and maintenance
- Affordable price
- Easy to clean
- For use with cura1 2719 Cordless Pad Transmitter
- Virtually no false alarms over the lifetime of the mat
- Made from heavy-duty non-slip material
- Cordless Floor Mat for the cordless range 2700

Buy this Cura1 2624 Cordless Floor Mat online using [NDIS Payment Method](#) or AfterPay Today.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au