

CURA1 2624 CORDLESS FLOOR MAT

\$375.00



SKU: 2624B



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Cura1 Wireless Floor Mat has the benefit of wireless out of bed detection methods in that there are no messy wires or potential trip hazards. As the resident places their foot on this floor mat the carer will be notified via the facilities nurse call system. This Cura1 brand mat has a “yoga” soft top and is waterproof and washable. It is low profile and made of a non-slip material. It has a press-stud protective pocket to house the transmitter/battery pack (sold separately).

Wireless floor mats and transmitter/battery packs are sold separately.

The transmitter/battery pack which is required for this mat is the [Cura1 Transmitter Module](#) and the Monitor is the [Cura1 Cordless Falls Monitor](#).

KEY FEATURES

- Simple operation and maintenance
- Affordable price
- Easy to clean
- For use with cura1 2719 Cordless Pad Transmitter
- Virtually no false alarms over the lifetime of the mat
- Made from heavy-duty non-slip material
- Cordless Floor Mat for the cordless range 2700

Buy this Cura1 2624 Cordless Floor Mat online using [NDIS Payment Method](#) or AfterPay Today.