



CURA1 CORDLESS PAD TRANSMITTER



SKU: CUR-2719

PRODUCT DESCRIPTION

Purpose & Benefits of the Cura1 Cordless Pad Transmitter

The Cura1 Cordless Pad Transmitter is designed to enhance patient monitoring systems by eliminating the need for cords between sensor pads or mats and monitoring devices. This wireless functionality reduces tripping hazards and simplifies setup, providing a safer and more efficient environment in healthcare settings.

Features

- **Seamless Integration:** Fits securely inside Cura1 cordless pads and mats, maintaining a discreet and unobtrusive profile.
- **Wireless Connectivity:** Operates at a 433 MHz frequency, offering a wireless range of 30 to 50 metres, depending on environmental conditions.
- **Long Battery Life:** Designed for low power consumption, providing up to one year of battery life under typical usage.
- **User-Friendly Indicators:** Incorporates status and in-use LED indicators, allowing caregivers to easily monitor the device's operational state.
- **Convenient Controls:** Features an on/off switch for straightforward operation and maintenance.

Specifications

- **Compatibility:** Compatible with Cura1 cordless pads and mats.
- **Operating Frequency:** 433 MHz.
- **Wireless Range:** Approximately 30 to 50 metres (subject to environmental factors).
- **Battery Life:** Up to one year under normal usage conditions.
- **Indicators:** LED status and in-use indicators.
- **Power Control:** On/off switch included.

The Cura1 Cordless Pad Transmitter offers a practical solution for enhancing patient safety and caregiver efficiency by providing reliable wireless communication between sensor pads and monitoring systems.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au