



## **CURA1 FALL PREVENTION BED/CHAIR PAD KIT**



**SKU:** CURA BED/CHAIR PAD KIT

# **PRODUCT DESCRIPTION**

The Cura1 Fall Prevention Bed/Chair Pad Kit is an all-inclusive solution designed to enhance safety for individuals at risk of falls. With both bed and chair monitoring capabilities, this kit offers comprehensive coverage to alert caregivers in various settings, ensuring quick responses and reducing the risk of accidents.

#### This kit includes:

- **Cura LED Pager:** Instantly receives alerts from the bed or chair pads, notifying caregivers of any movement or potential falls.
- Cura ActiveCare Pager Soft Cover: Protects the pager, offering durability and comfort for everyday use.
- SnapOn Bed Pad (760 x 250): Detects when an individual leaves the bed, alerting caregivers to prevent falls during night-time or rest periods.
- **SnapOn Pad Transmitter:** Wirelessly connects bed and chair pads to the pager, ensuring seamless, cordless communication.
- SnapOn Chair Pad (380 x 175): Triggers alerts when an individual leaves the chair, reducing the risk of falls or unassisted movement.

For further information go to **medimart.com.au** or give us a call: **1300 633 063** 



# SA

# **Brighton Showroom**

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

## **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

# **Reynella Showroom**

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

# **VIC**

#### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

# NT

#### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au

## **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au