



## CURA1 FALL PREVENTION BED/CHAIR PAD KIT

**\$686.00**



**SKU:** CURA BED/CHAIR PAD KIT

## PRODUCT DESCRIPTION

The Cura1 Fall Prevention Bed/Chair Pad Kit is an all-inclusive solution designed to enhance safety for individuals at risk of falls. With both bed and chair monitoring capabilities, this kit offers comprehensive coverage to alert caregivers in various settings, ensuring quick responses and reducing the risk of accidents.

### This kit includes:

- **Cura LED Pager:** Instantly receives alerts from the bed or chair pads, notifying caregivers of any movement or potential falls.
- **Cura ActiveCare Pager Soft Cover:** Protects the pager, offering durability and comfort for everyday use.
- **SnapOn Bed Pad (760 x 250):** Detects when an individual leaves the bed, alerting caregivers to prevent falls during night-time or rest periods.
- **SnapOn Pad Transmitter:** Wirelessly connects bed and chair pads to the pager, ensuring seamless, cordless communication.
- **SnapOn Chair Pad (380 x 175):** Triggers alerts when an individual leaves the chair, reducing the risk of falls or unassisted movement.