




## CURA1 LED PAGER 3508



**SKU:** CUR-3508

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Cura1 LED Pager 3508

The Cura1 LED Pager 3508 is a portable and versatile paging unit designed to enhance caregiver responsiveness by monitoring up to three Cura1 wireless devices simultaneously. Its compact design, combined with customisable alert settings, ensures that caregivers can promptly attend to patients' needs, thereby improving the quality of care provided.

### Features

- **Multi-Device Monitoring:** Capable of being programmed to monitor up to three wireless devices, allowing caregivers to manage multiple alerts efficiently.
- **Customisable Alerts:** Offers vibration and tone selection options, catering to caregivers who are hard of hearing or prefer silent alerts.
- **User-Friendly Interface:** Equipped with a cancel button for easy acknowledgment and management of alerts.
- **Portability:** Includes a belt clip for convenient carrying, ensuring caregivers can remain mobile while staying informed of any alerts.
- **Optional Accessories:** Compatible with an optional wall bracket for fixed positioning in care facilities.
- **Battery Operated:** Powered by two AA batteries (included), ensuring ease of maintenance and replacement.

### Specifications

- **Operating Frequency:** 433.92 MHz
- **Power Requirements:** 2 x AA batteries

- **Dimensions:** 83 mm (L) x 55 mm (W) x 26 mm (D)
- **Housing Material:** Durable construction suitable for daily use

The Cura1 LED Pager 3508 is an essential tool for caregivers, offering reliable monitoring and prompt alert notifications to ensure patient safety and well-being.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)