



CURA1 WIRELESS PAGER KIT

\$144.00



SKU: 4024

PRODUCT DESCRIPTION

Purpose & Benefits of the Cura1 Wireless Pager Kit

The Cura1 Wireless Pager Kit is a reliable and user-friendly communication system designed to enhance safety and responsiveness in care environments. Ideal for use in hospitals, aged care facilities, and homecare settings, this wireless kit allows caregivers to respond promptly to patient or resident calls. With its compact design and wireless functionality, it ensures seamless communication and peace of mind for both users and caregivers.

Features

- **Wireless Connectivity**: Offers reliable wireless communication between call buttons and the caregiver pager, eliminating the need for complex wiring.
- **Portable Caregiver Pager**: Compact and lightweight, the pager can be easily carried by caregivers, ensuring they remain responsive to calls at all times.
- **Multiple Call Buttons**: Supports up to 40 wireless devices, allowing multiple users to be connected to a single caregiver pager for efficient monitoring.
- **Audible and Vibration Alerts**: The pager features adjustable sound and vibration alerts, providing caregivers with flexible notification options to suit different environments.
- **Easy Installation**: The system is simple to set up with minimal technical expertise required, making it ideal for quick deployment in care facilities or homes.
- **Long Battery Life**: Designed with energy-efficient components, the pager, and call buttons offer extended battery life for dependable use.

Specifications

Transmitters use 1 x A23 12 volt battery (included)



- Pager uses 2 x AA 1.5 volt battery (included)
- Operating frequency 433 MHz
- Range adequate for most conventional home environments
- Durable ABS plastic construction
- 2 selectable alarm tones and off switch on the pager
- Green LED alarm indicator light

Cura1 Wireless Pager

[video width="1280" height="720" mp4="https://medimart.com.au/wp-content/uploads/2025/01/HC-Pager-Kit-Medimart.mp4"][/video]