


CUROCELL CIRRUS MATTRESS



SKU: MA0150-H



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The CuroCell Cirrus Mattress is an advanced therapeutic mattress designed for individuals requiring specialized care to prevent or treat pressure ulcers. This mattress is particularly suitable for use in hospitals, long-term care facilities, and home care settings. It combines innovative technology with user comfort to ensure optimal pressure redistribution and support.

FEATURES

- **Pressure Redistribution Technology:** The mattress utilizes a combination of air cells and foam to evenly distribute body weight. This technology helps to minimize pressure on vulnerable areas, significantly reducing the risk of pressure ulcer development.
- **Dynamic Air Support:** Equipped with an alternating pressure system, the CuroCell Cirrus dynamically adjusts the air pressure in the cells, providing regular movement and support changes that aid in pressure ulcer prevention and treatment.
- **Low Air Loss Feature:** The mattress incorporates a low air loss system that helps to manage the microclimate of the skin by allowing air to flow continuously, reducing moisture and enhancing skin integrity.
- **User-Friendly Controls:** Comes with an intuitive control unit that allows for easy adjustment of pressure settings, ensuring customized care tailored to the individual needs of each user.
- **Silent Pump Operation:** Designed to operate quietly, ensuring minimal disturbance to the user, which is especially beneficial in a home care setting or during nighttime use.
- **Comfort and Durability:** The top layer is made of high-quality, soft foam, providing additional comfort. The mattress materials are durable and designed for long-term use.

The CuroCell Cirrus Mattress is an excellent choice for patients with high-risk pressure ulcer profiles

or those in need of therapeutic support for existing skin integrity issues. Its advanced features and user-centric design make it a preferred option in advanced medical care settings, ensuring both comfort and effective treatment.

SPECIFICATIONS

- User weight up to 230kg
- Height: 21cm
- Warranty: 2 years

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au