

### **CUROCELL CIRRUS MATTRESS**



\$7,308.00 - \$8,431.00Price range: \$7,308.00 through \$8,431.00



**SKU:** MA0150-H

### PRODUCT DESCRIPTION

## **PURPOSE & BENEFITS**

The CuroCell Cirrus Mattress is an advanced therapeutic mattress designed for individuals requiring specialized care to prevent or treat pressure ulcers. This mattress is particularly suitable for use in hospitals, long-term care facilities, and home care settings. It combines innovative technology with user comfort to ensure optimal pressure redistribution and support.

### **FEATURES**

- **Pressure Redistribution Technology:** The mattress utilizes a combination of air cells and foam to evenly distribute body weight. This technology helps to minimize pressure on vulnerable areas, significantly reducing the risk of pressure ulcer development.
- **Dynamic Air Support:** Equipped with an alternating pressure system, the CuroCell Cirrus dynamically adjusts the air pressure in the cells, providing regular movement and support changes that aid in pressure ulcer prevention and treatment.
- Low Air Loss Feature: The mattress incorporates a low air loss system that helps to manage the microclimate of the skin by allowing air to flow continuously, reducing moisture and enhancing skin integrity.
- **User-Friendly Controls:** Comes with an intuitive control unit that allows for easy adjustment of pressure settings, ensuring customized care tailored to the individual needs of each user.
- **Silent Pump Operation:** Designed to operate quietly, ensuring minimal disturbance to the user, which is especially beneficial in a home care setting or during nighttime use.
- **Comfort and Durability:** The top layer is made of high-quality, soft foam, providing additional comfort. The mattress materials are durable and designed for long-term use.

The CuroCell Cirrus Mattress is an excellent choice for patients with high-risk pressure ulcer profiles



or those in need of therapeutic support for existing skin integrity issues. Its advanced features and user-centric design make it a preferred option in advanced medical care settings, ensuring both comfort and effective treatment.

# **SPECIFICATIONS**

• User weight up to 230kg

• Height: 21cm

• Warranty: 2 years