



THERAPEUTIC PILLOW COCCYX DIFFUSER CUSHION

~~\$120.00~~ Original price was: \$120.00. ~~\$105.00~~ Current price is: \$105.00. Sale Available Online Only! While Stocks Last



SKU: A14042001-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Coccyx Diffuser Cushion is a chair cushion designed to reduce coccyx pressure. The cut-out design alleviates tailbone pain, pressure and discomfort allowing the body to heal faster and more comfortably.

KEY FEATURES

- Orthopedically designed to reduce pressure on the coccyx or tailbone
- Cut-out suspends the coccyx weightlessly to relieve pressure
- Made from premium soft memory foam outer and heavy density traditional foam inner
- Heavy density traditional foam inner aids longevity and long term support
- Provides relief from bruising, pain or soreness
- 100% Australian made
- Ideal for wheelchairs or conventional seating

TWO COVER OPTIONS AVAILABLE

- Hospital grade SteriPlus polyurethane - moisture and bacteria resistant, features a convenient zip for easy removal, and can be wiped down with a cloth or removed and put in the washing machine.
- Soft faux-suede DuraFab cover - an attractive, soft, durable woven washable polyester fabric.

SPECIFICATIONS

Size

44 cm x 42 cm x 7.5 cm, 50 cm x 45 cm x 7.5 cm

Inner	100% open cell 'breathing' urethane foam and visco-elastic Memory Foam.
Dura-Fab Outer Cover	Woven polyester. It is a super-soft easy-clean washable fabric.
Steri-Plus Outer Cover	Moisture and bacteria-resistant polyurethane ply bonded to a knitted polyester backing.
Country of Manufacture	

Buy this Coccyx Wedge Cushion by TheraMed online at Medimart, with AfterPay available. Shop Now!