



ABLELIFE UNIVERSAL STAND ASSIST

\$237.00



SKU: 43001

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Universal Stand Assist makes it easy to stand from any couch or chair. The dual support handles make it easy to stand with maximum leverage. It adjusts to fit almost any couch or chair and blends in with your living room like another piece of furniture to maintain the comforting feeling of your home.

FEATURES

- Dual Support Handle makes sitting and standing easy
- · Adjustable -adjusts in height to fit most couches/chairs
- Easy installation snap together, slide it into the couch and adjust the height
- Low profile matches look and feel of your home

SPECIFICATIONS

- Handle Height Range from floor 54 80 cm
- Floor to Base Height 35 78cm
- Handle Height Adjustment from Base 19 26cm
- Width Adjustment 48 65cm
- Depth 60 68cm
- Weight Capacity 136kg
- Powder coated steel construction
- Weight of product 3kg