



DANMAR SOFT SHELL HELMET



SKU: DANSSHEINF

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This lightweight helmet is made from shock-absorbent foam that is fully ventilated. It's coated in vinyl, so it's easy to clean using just soap and water. With no hardware in the foam, there's no danger of the face connecting with anything hard or sharp. We feature a wide range of custom options, which we are happy to discuss with you so your individual needs are met. Comes with an easily adjustable snap-closure chinstrap.

FEATURES

- Helmet averages only 8 ounces
- Fully ventilated while providing head coverage
- Helmet is covered in vinyl, so it's easy to clean – only soap and water required!
- Comes with an easily adjustable snap-closure chinstrap.

SPECIFICATIONS

Since proper fit is important for the comfort of the person wearing headgear, measurements should be determined by professional personnel.

Measure the head circumference (A) at eyebrow level all the way around the head.

Occipital measurement (B) is taken from eyebrow level to the back of the head at the point that the helmet is to end.

Over the top of the Head (C) from the tip of one ear to the tip of the other ear.

SIZE CHART

(A) Circumference Range (all the way around)

	Soft Shell		Hard Shell	
	Inches	CM	Inches	CM
Infant	16.25 - 17	41.3 - 43.2	n/a	n/a
XX-Small	17.5 - 18.25	44.5 - 46.4	17.25 - 18	47.0 - 48.9
X-Small	19 - 19.75	48.3 - 50.2	19.25 - 20	48.9 - 50.8
Small	20.75 - 21.5	52.7 - 54.6	20.5 - 21.25	52.1 - 54.0
Medium	21.5 - 22.25	54.6 - 56.5	21.5 - 22.25	54.6 - 56.5
Large	23 - 23.75	58.4 - 60.3	22.5 - 23.25	57.2 - 59.1
X-Large	24 - 24.75	61 - 62.9	24.5 - 25.25	57.2 - 65.4
XX-Large	25.5 - 26.25	64.8 - 66.7	26 - 26.75	66.0 - 67.9

(B) Occipital

	Soft Shell		Hard Shell	
	Inches	CM	Inches	CM
Infant	11.5	29.2	n/a	n/a
XX-Small	12.75	32.4	12.75	32.4
X-Small	13.75	34.9	13.25	33.7
Small	14.5	36.8	15	38.1
Medium	15.25	38.7	16.5	41.9
Large	16.5	41.9	16.75	42.5
X-Large	17.5	44.5	17.5	44.5
XX-Large	18.5	47.0	18.75	47.6

(C) Ear To Ear

	Soft Shell		Hard Shell	
	Inches	CM	Inches	CM
Infant	8.25	21.0	n/a	n/a
XX-Small	9.5	24.1	9.25	23.5
X-Small	9.5	24.1	10	25.4
Small	10.5	26.7	11.75	29.8
Medium	11.5	27.9	12.5	31.8
Large	12	30.5	12.75	32.4
X-Large	13	33.0	13.25	33.7
XX-Large	13.75	34.9	15	38.1

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au