



DANMAR SOFT SHELL HELMET

\$412.00Online Only



SKU: DANSSHEINF

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This lightweight helmet is made from shock-absorbent foam that is fully ventilated. It's coated in vinyl, so it's easy to clean using just soap and water. With no hardware in the foam, there's no danger of the face connecting with anything hard or sharp. We feature a wide range of custom options, which we are happy to discuss with you so your individual needs are met. Comes with an easily adjustable snap-closure chinstrap.

FEATURES

- Helmet averages only 8 ounces
- Fully ventilated while providing head coverage
- Helmet is covered in vinyl, so it's easy to clean – only soap and water required!
- Comes with an easily adjustable snap-closure chinstrap.

SPECIFICATIONS

Since proper fit is important for the comfort of the person wearing headgear, measurements should be determined by professional personnel.

Measure the head circumference (A) at eyebrow level all the way around the head.

Occipital measurement (B) is taken from eyebrow level to the back of the head at the point that the helmet is to end.

Over the top of the Head (C) from the tip of one ear to the tip of the other ear.

SIZE CHART

(A) Circumference Range (all the way around)

	Soft Shell		Hard Shell	
	Inches	CM	Inches	CM
Infant	16.25 - 17	41.3 - 43.2	n/a	n/a
XX-Small	17.5 - 18.25	44.5 - 46.4	17.25 - 18	47.0 - 48.9
X-Small	19 - 19.75	48.3 - 50.2	19.25 - 20	48.9 - 50.8
Small	20.75 - 21.5	52.7 - 54.6	20.5 - 21.25	52.1 - 54.0
Medium	21.5 - 22.25	54.6 - 56.5	21.5 - 22.25	54.6 - 56.5
Large	23 - 23.75	58.4 - 60.3	22.5 - 23.25	57.2 - 59.1
X-Large	24 - 24.75	61 - 62.9	24.5 - 25.25	57.2 - 65.4
XX-Large	25.5 - 26.25	64.8 - 66.7	26 - 26.75	66.0 - 67.9

(B) Occipital

	Soft Shell		Hard Shell	
	Inches	CM	Inches	CM
Infant	11.5	29.2	n/a	n/a
XX-Small	12.75	32.4	12.75	32.4
X-Small	13.75	34.9	13.25	33.7
Small	14.5	36.8	15	38.1
Medium	15.25	38.7	16.5	41.9
Large	16.5	41.9	16.75	42.5
X-Large	17.5	44.5	17.5	44.5
XX-Large	18.5	47.0	18.75	47.6

(C) Ear To Ear

	Soft Shell		Hard Shell	
	Inches	CM	Inches	CM
Infant	8.25	21.0	n/a	n/a
XX-Small	9.5	24.1	9.25	23.5
X-Small	9.5	24.1	10	25.4
Small	10.5	26.7	11.75	29.8
Medium	11.5	27.9	12.5	31.8
Large	12	30.5	12.75	32.4
X-Large	13	33.0	13.25	33.7
XX-Large	13.75	34.9	15	38.1