



DAYS 102 ROLLATOR

\$190.00 Original price was: \$190.00.<u>\$171.00</u>Current price is: \$171.00.



SKU: DAYS-102-RED

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Days 102 Series, where lightweight meets strength. Crafted with a high-quality aluminium frame, this design offers the perfect balance between a lightweight structure and enduring build quality.

FEATURES

- **Durable 6" Wheels:** Navigate smoothly on various terrains.
- Ergonomic Handgrips: Designed for comfortable and easy gripping.
- Adjustable Handle Height: Customise your grip with 5 different settings at intervals of 25mm.
- **Rest Seat:** Take a moment to sit and relax during your outings.
- Storage Bag: Conveniently store your essentials on the go.
- **Secure Braking System:** Featuring loop lockable brakes for enhanced security and convenience. Plus, adjust the brake tension effortlessly to fit your needs.
- Foldable Backrest: Ensures comfort and easy storage.
- Easy Fold Mechanism: No hassle storage and transport. Simply pull the seat strap and fold it.
- Lightweight: Weighing only 6kg, it's easy to manoeuvre and transport.
- Strong Weight Limit: Supports a safe working load (SWL) of up to 125kg.

SPECIFICATIONS

Rollator Type4 WheelFrame MaterialAluminiumSeat Height Front500 mm



500 mm

Seat Height Rear Seat Width

Seat Width350 mmSeat Depth320 mmOverall Depth640 mm

Overall Height 755-875 mm

Gross Weight 6 kg **Max SWL** 125 kg

Brand Days Healthcare