



## **DAYS 111 BARIATRIC ROLLATOR**



**SKU:** PAT-7500348

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The Days 111 Bariatric Rollator is a high-strength mobility aid designed to support larger individuals with stability, comfort, and confidence. Constructed from reinforced steel, it helps users maintain balance, reduce fall risk, and move independently across both indoor and outdoor surfaces.

This rollator features a padded backrest and secure seat for comfortable resting, a built-in safety bar to prevent tipping backward while seated, and a front storage basket for carrying personal items. It folds down easily for transport and storage, making it a reliable solution for home, healthcare, or rehabilitation use.

### **Features**

- Heavy-duty steel frame supports up to 317.5 kg
- Padded seat and backrest for comfortable resting
- Rear safety bar prevents tipping when seated
- Height-adjustable handles with loop-style brakes
- 8" wheels for smooth navigation indoors and outdoors

- Large under-seat basket for storage
- Foldable for convenient transport and storage
- Ideal for bariatric users at home or in care facilities

## Specifications

- **Maximum User Weight:** 317.5 kg
- **Seat Width:** 49 cm
- **Seat Depth:** 35 cm
- **Seat Height:** 55 cm
- **Handle Height Range:** 83-96.5 cm
- **Overall Width:** 61 cm
- **Overall Depth:** 71 cm
- **Product Weight:** 13 kg
- **Wheel Size:** 8"
- **Frame Material:** Heavy-duty steel
- **Colour:** Red

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

**SA**

**Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

**Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

**Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

**VIC**

**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT**

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)