



## **DAYS 111 BARIATRIC ROLLATOR**

**\$460.00**



**SKU:** PAT-7500348

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The Days 111 Bariatric Rollator is a high-strength mobility aid designed to support larger individuals with stability, comfort, and confidence. Constructed from reinforced steel, it helps users maintain balance, reduce fall risk, and move independently across both indoor and outdoor surfaces.

This rollator features a padded backrest and secure seat for comfortable resting, a built-in safety bar to prevent tipping backward while seated, and a front storage basket for carrying personal items. It folds down easily for transport and storage, making it a reliable solution for home, healthcare, or rehabilitation use.

### **Features**

- Heavy-duty steel frame supports up to 317.5 kg
- Padded seat and backrest for comfortable resting
- Rear safety bar prevents tipping when seated
- Height-adjustable handles with loop-style brakes
- 8" wheels for smooth navigation indoors and outdoors

- Large under-seat basket for storage
- Foldable for convenient transport and storage
- Ideal for bariatric users at home or in care facilities

## Specifications

- **Maximum User Weight:** 317.5 kg
- **Seat Width:** 49 cm
- **Seat Depth:** 35 cm
- **Seat Height:** 55 cm
- **Handle Height Range:** 83-96.5 cm
- **Overall Width:** 61 cm
- **Overall Depth:** 71 cm
- **Product Weight:** 13 kg
- **Wheel Size:** 8"
- **Frame Material:** Heavy-duty steel
- **Colour:** Red