



DAYS ALUMINIUM FOLDING SHOWER CHAIR

\$142.00



SKU: JAN-BE04

PRODUCT DESCRIPTION

Purpose & Benefits of the Days Aluminium Folding Shower Chair

The Days Aluminium Folding Shower Chair is designed to provide a safe and comfortable seating option for individuals who need extra support while showering. Its lightweight, foldable design makes it easy to store and transport, making it ideal for use in both residential and care settings. The chair's sturdy aluminium frame is rust-resistant, ensuring durability and longevity in wet environments. It offers enhanced safety and comfort, allowing users to maintain independence and enjoy a more secure showering experience.

Key Features

- **Lightweight, Foldable Design:** The chair folds flat for easy storage and portability, making it convenient for travel or when space is limited.
- **Sturdy Aluminium Frame:** Made from rust-resistant aluminium, the chair provides reliable support and is suitable for long-term use in wet environments.
- **Height Adjustable Legs:** The chair's legs can be adjusted to accommodate different user heights, ensuring a comfortable and customised seating position.
- **Non-Slip Rubber Feet:** Equipped with non-slip rubber feet for enhanced stability, the chair prevents slipping and ensures safety during use.
- **Comfortable Plastic Seat:** Features a contoured plastic seat with drainage holes, providing comfort while allowing water to flow through, preventing water build-up.

Specifications

- **Seat Width:** 41 cm

- **Height Range:** Adjustable from 47 cm to 57 cm
- **Weight Capacity:** 130 kg
- **Material:** Rust-resistant aluminium frame with a durable plastic seat
- **Cleaning:** Wipe clean with a damp cloth or mild detergent