



## DAYS ALUMINIUM FOLDING SHOWER CHAIR



**SKU:** JAN-BE04

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Days Aluminium Folding Shower Chair

The Days Aluminium Folding Shower Chair is designed to provide a safe and comfortable seating option for individuals who need extra support while showering. Its lightweight, foldable design makes it easy to store and transport, making it ideal for use in both residential and care settings. The chair's sturdy aluminium frame is rust-resistant, ensuring durability and longevity in wet environments. It offers enhanced safety and comfort, allowing users to maintain independence and enjoy a more secure showering experience.

### Key Features

- **Lightweight, Foldable Design:** The chair folds flat for easy storage and portability, making it convenient for travel or when space is limited.
- **Sturdy Aluminium Frame:** Made from rust-resistant aluminium, the chair provides reliable support and is suitable for long-term use in wet environments.
- **Height Adjustable Legs:** The chair's legs can be adjusted to accommodate different user heights, ensuring a comfortable and customised seating position.
- **Non-Slip Rubber Feet:** Equipped with non-slip rubber feet for enhanced stability, the chair prevents slipping and ensures safety during use.
- **Comfortable Plastic Seat:** Features a contoured plastic seat with drainage holes, providing comfort while allowing water to flow through, preventing water build-up.

### Specifications

- **Seat Width:** 41 cm

- **Height Range:** Adjustable from 47 cm to 57 cm
- **Weight Capacity:** 130 kg
- **Material:** Rust-resistant aluminium frame with a durable plastic seat
- **Cleaning:** Wipe clean with a damp cloth or mild detergent

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)