




DAYS THE MACK BARIATRIC WALKER

\$418.00 - \$449.00



SKU: NOV-

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Days Mack Bariatric Walker is constructed from heavy-duty aluminium for durability and features a removable padded backrest with easy release pins. Handles are ergonomically designed and the handbrakes can be locked for stability. The padded seat allows for comfort and the shopping basket is a handy aid and can be removed when folding for transportation or storage.

SPECIFICATIONS

Supa Mack

- Weight capacity: 225 kg
- Product weight: 9.4 kg
- Castor wheel front/rear: 8" x 2" solid
- Seat width: 610 mm
- Seat floor height: 550 mm

Name	Low Mack	High Mack
Overall Width	700mm	700mm
Overall Depth	640mm	640mm
Floor to Seat Height	490mm	580mm
Seat Width	500mm	500mm

Handle Height Range	750mm - 860mm	870mm - 980mm
Weight of Walker	9.4kg	9.4kg
Colour	Red	
Safe Working Load	150kg	