



DAYS THE MACK BARIATRIC WALKER

\$400.00



SKU: NOV-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Days Mack Bariatric Walker is constructed from heavy-duty aluminium for durability and features a removable padded backrest with easy release pins. Handles are ergonomically designed and the handbrakes can be locked for stability. The padded seat allows for comfort and the shopping basket is a handy aid and can be removed when folding for transportation or storage.

SPECIFICATIONS

Supa Mack

Weight capacity: 225 kgProduct weight: 9.4 kg

Castor wheel front/rear: 8" x 2" solid

• Seat width: 610 mm

• Seat floor height: 550 mm

Name	Low Mack	High Mack
Overall Width	700mm	700mm
Overall Depth	640mm	640mm
Floor to Seat Height	490mm	580mm
Seat Width	500mm	500mm
Handle Height Range	750mm - 860mm	870mm - 980mm



Weight of Walker 9.4kg 9.4kg

Colour Red

Safe Working Load 150kg