



## DAYS THE MACK BARIATRIC WALKER

**\$400.00 - \$449.00** Price range: \$400.00 through \$449.00



**SKU:** NOV-

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Days Mack Bariatric Walker is constructed from heavy-duty aluminium for durability and features a removable padded backrest with easy release pins. Handles are ergonomically designed and the handbrakes can be locked for stability. The padded seat allows for comfort and the shopping basket is a handy aid and can be removed when folding for transportation or storage.

### SPECIFICATIONS

#### Supa Mack

- Weight capacity: 225 kg
- Product weight: 9.4 kg
- Castor wheel front/rear: 8" x 2" solid
- Seat width: 610 mm
- Seat floor height: 550 mm

Name	Low Mack	High Mack
<b>Overall Width</b>	700mm	700mm
<b>Overall Depth</b>	640mm	640mm
<b>Floor to Seat Height</b>	490mm	580mm
<b>Seat Width</b>	500mm	500mm
<b>Handle Height Range</b>	750mm - 860mm	870mm - 980mm

<b>Weight of Walker</b>	9.4kg	9.4kg
<b>Colour</b>	Red	
<b>Safe Working Load</b>	150kg	