



# DAYS BARIATRIC WALKER WITH DOUBLE 5" WHEELS

\$285.00



**SKU:** PAT-081561653

### PRODUCT DESCRIPTION

## **Purpose & Benefits**

The **Days Bariatric Walker with Double 5" Wheels** provides strong, stable, and lightweight mobility support for bariatric users. Constructed from durable aluminium with double-cross bracing, it offers excellent stability without unnecessary bulk, while its dual-release folding mechanism ensures safe and easy operation.

#### **Features**

- Lightweight aluminium frame designed for maximum strength and support
- Dual-release buttons provide audible locking cues for secure adjustments and folding
- Double-cross bracing enhances frame stability and durability
- Contoured hand grips ensure comfortable use and reduced strain
- Fitted with double 5" front wheels for smooth mobility
- Rubber-like polyurethane tips provide secure traction and stability
- Folds easily for convenient transport and storage



# **Specifications**

• Material: Aluminium

• Wheel Size: Double 5"

• Height Adjustment Range: Suitable for users 1.6 m - 1.85 m

• Overall Width: 28" (71 cm)

• Width Between Hand Grips: 20" (51 cm)

• Width Inside Base of Legs: 24" (61 cm)

• **Depth:** 16" (41 cm)

• Product Weight: Approx. 4 kg

• Maximum User Weight: 317.5 kg