




## **DAYS BARIATRIC WALKER WITH DOUBLE 5" WHEELS**



**SKU:** PAT-081561653

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The **Days Bariatric Walker with Double 5" Wheels** provides strong, stable, and lightweight mobility support for bariatric users. Constructed from durable aluminium with double-cross bracing, it offers excellent stability without unnecessary bulk, while its dual-release folding mechanism ensures safe and easy operation.

### **Features**

- Lightweight aluminium frame designed for maximum strength and support
- Dual-release buttons provide audible locking cues for secure adjustments and folding
- Double-cross bracing enhances frame stability and durability
- Contoured hand grips ensure comfortable use and reduced strain
- Fitted with double 5" front wheels for smooth mobility
- Rubber-like polyurethane tips provide secure traction and stability
- Folds easily for convenient transport and storage

## Specifications

- **Material:** Aluminium
- **Wheel Size:** Double 5"
- **Height Adjustment Range:** Suitable for users 1.6 m – 1.85 m
- **Overall Width:** 28" (71 cm)
- **Width Between Hand Grips:** 20" (51 cm)
- **Width Inside Base of Legs:** 24" (61 cm)
- **Depth:** 16" (41 cm)
- **Product Weight:** Approx. 4 kg
- **Maximum User Weight:** 317.5 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)