




DAYS COMFORT CHAIR



SKU: ANI-DAYCOMCA1

PRODUCT DESCRIPTION

Purpose & Benefits of the Days Comfort Chair

The Days Comfort Chair is designed to provide exceptional support and comfort for individuals requiring extended periods of seated care. Ideal for elderly users or those with specific postural and pressure care needs, this chair offers a fully adjustable design to enhance both user comfort and caregiver convenience.

Features

- **Adjustable Recline** – Offers a recline range from 90° to 170°, allowing users to find their optimal seating or resting position.
- **Customisable Support** – Features an adjustable backrest, calf rest, neck support, and footplates to suit individual requirements.
- **Pressure Relief System** – Integrated cushioned support with adjustable air cells to aid in pressure care management.
- **Facilitated Transfers** – All four lateral support wings fold away, making transfers easier for caregivers.
- **Integrated Dining Tray** – Includes a removable tray for meal times, activities, or reading.
- **Durable Construction** – Built with a stainless steel frame and covered in breathable, fire-retardant polyurethane material.
- **Enhanced Mobility** – Equipped with three lockable castors and one steering castor for secure positioning and smooth manoeuvrability.
- **Convenient Storage** – Includes a built-in magazine pocket and footplates for added practicality.

Specifications

- **Total Length:** 350 mm – 1880 mm
- **Total Width:** 730 mm
- **Seat Height:** 630 mm
- **Seat Width:** 490 mm
- **Seat Depth:** 500 mm
- **Back Recline:** 102° – 136°
- **Seat Tilt:** 0° – 18°
- **Leg Elevation:** Down: 70° | Up: 5°
- **Table Size:** 400 mm x 405 mm x 20 mm
- **Table Slide:** 230 mm
- **Table to Floor Height:** 880 mm
- **Product Weight:** 43.5 kg
- **Safe Working Load:** 120 kg

The Days Comfort Chair is a reliable and versatile seating solution, offering superior comfort and postural support for individuals in home care and healthcare environments.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au