



DAYS GUTTER WALKER WITH HANDBRAKES



SKU: DAYS-203B

PRODUCT DESCRIPTION

Purpose & Benefits of the Days Gutter Walker

The Days Gutter Walker is specially designed to provide support and mobility assistance to individuals who need forearm support due to reduced hand and wrist strength or those who cannot use standard walkers due to arthritis, neurological disorders, or other conditions affecting their ability to grip. This walker features padded gutter-style forearm supports that allow the user to lean on their forearms instead of relying on hand grip strength, promoting an upright posture and a more natural walking pattern.

Features

- **Forearm Support:** Equipped with padded gutters to comfortably support the user's forearms, reducing the strain on hands and wrists and enabling easier mobility.
- **Adjustable Height:** The height of both the handles and the forearm gutters can be adjusted to fit the user's specific needs, ensuring comfort and optimal support.
- **Sturdy Frame:** Constructed from high-quality materials that ensure durability and stability, supporting safe and secure movement.
- **Easy Maneuverability:** Designed with wheels that facilitate smooth movement across various surfaces, ideal for indoor and outdoor use.
- **Foldable Design:** Can be easily folded, making it convenient for storage and transportation.
- **Brake System:** Features an effective brake system that provides additional safety by allowing the user to stop the walker easily when needed.

Specifications

- **Safe Working Load:** 150kg

- **Dimensions:**

- Overall Height: 1100 - 1340mm
- Overall Width: 630mm
- Overall Depth: 900mm
- Overall Forearm Height: 980 - 1220mm
- Width Between Wheels: 580mm
- Wheels: 190mm (7.48")

- **Care Instructions:** Easy to maintain; regular cleaning with a damp cloth and routine checks of the braking system and structural integrity are recommended.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au